

EYEPIECE GRATICULES (RETICLES)


METALLURGY

Standard pattern discs for metallurgical stereometric analysis of grain size in polished metal sections.

ASTM Austenite 1:1 Grain Sizing Disc – G41


Reference: VDEH 1510-61

Pattern	Description	Diameter	Order Code
G41	ASTM Grain sizing austenite. Glass sandwich.	19mm	01A19064
		21mm	01A21064
		Special	01ASP064




ASTM E112 Plate 1 Grain Sizing Disc – G42

Pattern	Description	Diameter	Order Code
G42	Grain sizing E112. Glass sandwich.	19mm	01A19065
		21mm	01A21065
		Special	01ASP065



ASTM Carbide Grain Sizing Chart – G43


Pattern	Description	Diameter	Order Code
G43	ASTM Grain sizing carbide. Glass sandwich.	19mm	01A19066
		21mm	01A21066
		Special	01ASP066



ASTM E45 - G44

Reference: ASTM E45

Pattern	Description	Diameter	Order Code
G44	ASTM E45. Root 2 sides. 7.1mm square. 10mm scale. Glass sandwich.	19mm	01A19086
		21mm	01A21086
		Special	01ASP086




We manufacture components (encoder discs, graticules, metal foils, resolution standards) to custom design, if it is not in the catalogue please contact us for pricing information.

EYEPIECE GRATICULES (RETICLES)

METALLURGY


ASTM E19-46 Grain Sizing Disc – G45

Pattern	Description	Diameter	Order Code
G45	ASTM Grain sizing E19-46. Glass sandwich.	19mm 21mm Special	01A19067 01A21067 01ASP067



ASTM E19-46 Grain Sizing Disc Root 2 – G46


Pattern	Description	Diameter	Order Code
G46	ASTM Grain sizing E19-46 Root 2. Glass sandwich.	19mm 21mm Special	01A19068 01A21068 01ASP068



Circular Grid ASTM 24 points – G54

Reference: ASTM E562

Pattern	Description	Diameter	Order Code
G54	ASTM 24 point circular grid. Glass sandwich.	19mm 21mm Special	01A19078 01A21078 01ASP078



Square Grid 25 points – G55

Reference: ASTM E562.

Pattern	Description	Diameter	Order Code
G55	ASTM 25 point. Square grid. Glass sandwich.	19mm 21mm Special	01A19079 01A21079 01ASP079

