


Mini-Fiber Optics Holders

- Designed to separately secure clad and jacketed fiber using V-grooves and magnet.
- There are three diameter versions, 250, 500 and 900µm to hold standard jacketed fiber sizes
- The V-groove for the clad fiber will secure diameters from 60 to 250µm
- The holders are designed to be used in combination with the Mini-Fiber adapters



MINI-FIBER OPTICS HOLDER

Mini Fiber Holders

Clad Diameter ØA (µm)	Jacket Diameter ØB (µm)	Price	PART NUMBER
60 ~ 130	150 ~ 250		111-0701
125 ~ 250	500.0		111-0702
125 ~ 250	900.0		111-0703

ORDERING
&
TECHNICAL SUPPORT
(949) 851-5881
FAX (949) 851-5058
E-MAIL
sales@optosigma.com
WEB
www.optosigma.com




Adapters for Mini-Fiber Optics


- Designed for use mini optical fiber holders
- Used to secure a small optical fiber holder to a stage or holder
- Offered in three versions to suit your application




Part No. 111-0673



Part No. 111-0672



Part No. 111-0671



ORDERING
&
TECHNICAL SUPPORT

(949) 851-5881
FAX (949) 851-5058

E-MAIL
sales@optosigma.com

WEB
www.optosigma.com

Adapter for Mini-Fiber Optics Holders

Description	Price	PART NUMBER
for mounting flat surface (M2, M3 counterbored)		111-0671
for mounting flat surface (M2, M3 counterbored)		111-0672
for mounting flat surface (M2, M3 counterbored)		111-0673

